



Purpose: Basic checklist for Krios operators to set up sessions for users

Authors:

Approved:

Date: 2020/07/08

1. Purpose:

- 1.1. To setup the Krios system and Legimon ahead of users starting Legimon Remote.

2. Scope:

- 2.1. Checklist that should be reviewed by microscope operator ahead of giving control of data collection to a Legimon Remote users.

3. Definitions

- 3.1. Legimon is a system designed for automated collection of images from a transmission electron microscope; it includes the python-side programs written in python and c, the MySQL database and server, and the mainly php-based image and data viewers on a web server.
- 3.2. Legimon Remote is an infallible web interface that allows center staff operators to hand over control of the requested transmission electron microscope to users for their own data acquisition.

4. Responsibilities:

4.1. Start of Day:

- 4.1.1. Setup a Legimon session
- 4.1.2. Determine eucentric height/focus
- 4.1.3. Find vacuum and tune energy filter
- 4.1.4. After an annealing cycle: GMS gains
- 4.1.5. Hardware dark & Legimon gain reference
- 4.1.6. New dose for exposure
- 4.1.7. Set up ice thickness node
- 4.1.8. Beam tilt pivot point and rotation centering
- 4.1.9. Objective aperture centering and astigmatism correction
- 4.1.10. Coma correction
- 4.1.11. Preset alignment
- 4.1.12. Correct image shift coma effect=true
- 4.1.13. Test images
- 4.1.14. Repeat if necessary
- 4.1.15. Switch to user's grid and collect an atlas

4.2. End of Day:

- 4.2.1. Toggle error notification ON
- 4.2.2. Check enn/esn/fa/fc beam centered over camera
- 4.2.3. Objective aperture is inserted (unless no requested)
- 4.2.4. Random defocus range is set
- 4.2.5. Final queue is submitted
- 4.2.6. Frame alignment and CTF estimation running -a/-a-DW
- 4.2.7. Screensaver and screen lock is off
- 4.2.8. Legimon station is VNC capable
- 4.2.9. Previous user has signed off, i.e. no "switch user"
- 4.2.10. Can remote desktop to microscope, camera and Legimon workstation
- 4.2.11. Data transfer started Globus/external drive
- 4.2.12. Remind user of grid storage/pick up