

CRYOEM 001 : CRYO SAMPLE PREPARATION

NCCAT Embedded Training — Master Class series

June 16-17, 2021

NATIONAL CENTER FOR
CRYOEM ACCESS & TRAINING



New York Structural
Biology Center

SIMONS ELECTRON
MICROSCOPY CENTER



WEEK 2

JUNE 6/16-17

Wed 6/16 3-5pm



Overview of cryo sample prep



Introduction to merit badges



Practical: Vitrobot Mark IV and Glacios

Thur 6/17 2-5pm



Screening microscope (Glacios) in detail



Screening with Leginon

Core
knowledge

Sample
preparation

Microscope
operations &
Data
collection

Processing
&
Data analysis

SCHEDULE

I. Sample purification and grid preparation

a) cryoEM merit badges

b) Chameleon demo

overview of blot free vitrification vs plunge freezing methods

II. Grid screening & evaluation

a) Sample holders

-Side entry systems: Gatan 626/Elsa holder and loading

-Autoloader systems: autogrid clipping and loading

b) F20 setup and demo of screening with Leginon

III. Cryo-EM data collection

a) Glacios setup and advanced sample

screening/preliminary data collection with Leginon

b) Krios high res data collection with Leginon

IV. Image (pre)-processing

a) On the fly feedback cryoSPARC live

b) Working with your own data

NCCAT CROSS-TRAINING RESOURCES

Supplemental cryoEM masterclass materials for cohort3

- Intro to sample preparation – negative stain and screening on an F20
- Negative stain and F20 Manuals and Checklists:
 - Solarus plasma cleaner manual
 - Negative-stain manual
 - Negative-Stain Independence Checklist
 - TF20 user manual
 - T12-F20 Checklist



CRYOEM MERIT BADGES



Broadening access to high-resolution
cryo-electron microscopy and tomography

HOME CryoEM CENTERS CURRICULUM DEVELOPMENT CryoET CENTERS EVENTS RESOURCES

[< return to all badges](#)

TFS Vitrobot Mark IV

Category: Sample preparation
Sub-category: Cryogenic work

Plunge freezing and instrument certification for Vitrobot Mark IV.

+ Essential base knowledge

+ Knowledge quiz

+ Center Specific Policies

+ Demonstration

+ Supervised Practice

+ Practical test




Figure 1. Vitrobot assembled and turned on. **A)** Screen. **B)** Environmental chamber with blotting pads. **C)** Humidifier. **D)** ethane lift.




Figure 2. Ethane Holder. **A)** Spider. **B)** Brass Ethane Cup. **C)** Gridbox Holder. **D)** Base / Liquid Nitrogen Container. **E)** Anti-contamination Ring.

Recertification period

- Sample preparation merit badges are valid for ~1yr.
- Recertification (to maintain active status) requires passing the practical test with one center staff member. If supervised training is needed to pass the practical test, this can be arranged.

CRYOEM MERIT BADGES

<https://cryoem101.org/selftest/?test=19>



CryoEM
TRANSFORMATIVE HIGH RESOLUTION
CRYO-ELECTRON MICROSCOPY

Broadening access to high-resolution
cryo-electron microscopy and tomography

[HOME](#) [CryoEM CENTERS](#) [CURRICULUM DEVELOPMENT](#) [CryoET CENTERS](#) [EVENTS](#) [RESOURCES](#)

[< return to all badges](#)

TFS Vitrobot Mark IV

Category: Sample preparation
Sub-category: Cryogenic work

Plunge freezing and instrument certification for Vitrobot Mark IV.

+ Essential base knowledge

+ Knowledge quiz

+ Center Specific Policies

+ Demonstration


+ Supervised Practice

+ Practical test



Recertification period

- Sample preparation merit badges are valid for ~1yr.
- Recertification (to maintain active status) requires passing the practical test with one center staff member. If supervised training is needed to pass the practical test, this can be arranged.



Begin Quiz: Merit Badge Knowledge Quiz – TFS Vitrobot Mark IV

When you're ready, fill in your information and click the "Start the Quiz" button

Test of foundational knowledge for Vitrobot use. You must answer 20 of the 23 questions correctly to pass. You may take the quiz multiple times.

First Name

Last Name




E-Mail Address

Affiliation

Position


[Start the Quiz](#)

PLUNGE FREEZING




PLUNGE FREEZING


- Liquid ethane is a suitable coolant.
- Liquid nitrogen boils on contact, which makes it a poor coolant for cryo-EM.
- Cooling speed faster than 10^5 - 10^6 K/s ensure the formation of vitrified ice.



Setup of liquid ethane
(Image from Wen Jiang)





Cooling speed &
forms of ice




Different forms of ice contamination




WHAT DO GRIDS LOOK LIKE?




WHAT ISSUES ARISE?




Aggregating in ice




Preferred orientation



Particles not going into holes



Rejecting 90% of particles



Particles disappearing in ice


WHAT ISSUES ARISE?




Hemagglutinin



Hemagglutinin



GDH



T20S Proteasome


110 nm
ice



Aldolase

45 nm
ice


125 nm
ice



Aldolase

50 nm
ice


180 nm
ice



Apoferritin

135 nm
ice

115 nm
ice



DNAB Helices

110 nm
ice

Noble AJ, et al.
Routine single
particle CryoEM
sample and grid
characterization
by tomography.
Elife. 2018;7.



Alex Noble

REAGENTS FOR IMPROVING VITRIFICATION OF CRYO-EM GRIDS USED IN SINGLE PARTICLE ANALYSIS.

Surfactants and Cryoprotectants	Amount	Conc.	CMC	Class
Fluorinated Octyl Maltoside (FOM)	100 µl	0.41% (w/v)	0.07% (w/v)	non-ionic detergent
Hexadecyl-trimethyl-ammonium Bromide (CTAB)	100 µl	0.34% (w/v)	0.03% (w/v)	cationic detergent
n-Decyl-β-D-Maltoside (DM)	100 µl	0.87% (w/v)	0.09% (w/v)	non-ionic detergent
n-Decyl-α-D-Maltoside (DaM)	100 µl	0.46% (w/v)	0.08% (w/v)	non-ionic detergent
n-Dodecyl-β-D-Maltoside (DDM)	100 µl	0.09% (w/v)	0.01% (w/v)	non-ionic detergent
Sodium Deoxycholate	100 µl	1.66% (w/v)	0.17% (w/v)	anionic detergent
Triton X-100	100 µl	0.15% (w/v)	0.01% (w/v)	non-ionic detergent
Tween 20	100 µl	1% (w/v)	0.01% (w/v)	non-ionic detergent
CHAPSO	100 µl	2.5% (w/v)	0.5% (w/v)	zwitterionic detergent
Amphipol A8-35	100 µl	5% (w/v)		anionic surfactant
Glycerol	1 ml	30% (w/v)		cryoprotectant

- [1] Noble *et al.* (2018) Routine Single Particle CryoEM Sample and Grid Characterization by Tomography. DOI: 10.7554/eLife.34257.
- [2] Thonghin *et al.* (2018) Cryo-electron microscopy of membrane proteins. *Methods* **147**:176.
- [3] Drulyte *et al.* (2018) Approaches to altering particle distributions in cryo-electron microscopy sample preparation. *Acta Cryst. D* **74**:560.
- [4] Glaeser *et al.* (2017) Opinion: hazards faced by macromolecules when confined to thin aqueous films. *Biophys Rep* **3**:1.
- [5] Gatsogiannis *et al.* (2016). Membrane insertion of a Tc toxin in near-atomic detail. *Nat. Struct. Mol. Biol.* **23**:884.
- [6] Efremov *et al.* (2015) Architecture and conformational switch mechanism of the ryanodine receptor. *Nature* **517**:39.

<https://www.mitegen.com/product/cryo-em-vitrification-starter-kit/>

REAGENTS FOR IMPROVING VITRIFICATION OF CRYO-EM GRIDS USED IN SINGLE PARTICLE ANALYSIS.

PDB Release Date	PDB	Protein	Additive
2020-01-08	6PWN	MscS mechanosensitive channel	0.01% f-OM
2019-09-04	6KG7	Piezo2 mechanosensitive channel	0.65 mM f-FC8
2019-08-28	6QTI	Nicotinamide nucleotide proton channel	0.05% CHAPS
2019-08-07	6R7L	SecYEG translocon	0.2% f-OM
2019-02-06	6E0H	TMEM16 scramblase	3 mM f-FC8
2018-12-19	6N3Q	Sec protein-translocation channel complex	3 mM f-FC8
2018-11-07	6H3I	Type 9 secretion system translocon	1.5 mM f-FC8 or 0.7 mM f-OM
2018-10-24	6DMR	TRPV5 ion channel	3 mM f-FC8
2018-10-17	6D3R	CFTR	3 mM f-FC8
2018-09-26	6HJR	Influenza Hemagglutinin	2% Octyl Glucoside
2018-08-08	6FOO	Ryanodine receptor 1	0.2% f-OM
2018-08-01	6CJQ	SthK CNG Potassium channel	3 mM f-FC8
2018-05-23	5YX9	TRPC6 ion channel	0.5 mM f-OM
2018-01-31	6C0V	P-Glycoprotein transporter ABCB1	3 mM f-FC8
2017-12-27	6B5V	TRPV5 ion channel	3 mM f-FC8
2017-12-13	6BPQ	TRPM8 channel	2% DMSO

Glaeser, RM, et al. (2017) Biophys Rep 3(1), 1-7.


Noble, AJ, et al. (2018) Nat Methods 15(10), 793-795.

Drulyte, I et al. (2018) Acta Crystallogr D Struct Biol 74(Pt 6), 560-571.

Chen, J, et al. (2019) J Struct Biol X Volume 1. DOI: 10.1016/j.yjsbx.2019.100005

<https://www.anatrace.com/Landing/2020/Mar20-Newsletter>

DIFFICULT SPECIMENS



Small protein

- VPP
- Thinner ice

Protein denaturation/Dissociation of protein complex

Continuous carbon film

Graphene oxide

Cross-linking (GraFix)

Preferred orientation

Tilt stage

Cross-linking

Detergent

Glow-discharging conditions

Support film (Graphene oxide)

Image analysis (3D classification)

Flexibility

Focused classification (subtraction)

Multibody refinement

Filamentous protein


- Segmented analysis

VI. Low concentration

Multiple blots

Affinity grids

EMGWEB.NYSBC.ORG Appion and Leginon Tools



EMGWEB.NYSBC.ORG

ninfo / nccat032020

002

CRYOEM 101

START
HERE..

WEEK 1
JUNE 6/09-10

WEEK 2
JUNE 6/16-17

WEEK 3
WK OF JUNE 21+



Practical: High end data collection



Data collection on Krios



Leginon automation



Pre-processing / cryoSPARC LIVE